



Atomic Force Microscopy's Characterization of Graphene Nanoplate and bubble




Jinyu Chong[†], Yuxuan Xue[§], Ning Xi^{§, *}

[†] Department of Chemistry, Soochow University, Suzhou 215123, P.R. China

[§] Department of Industrial and Manufacturing Systems Engineering, The University of Hong Kong, Hong Kong, China

Introduction

- Atomic Force Microscopy (AFM) has been widely used since it was first invented in 1980s. Through a sharp tip scanned back and forth across a sample's surface, AFM determine the surface topology of nanomaterials with high resolution. There is also a static mode to analyze force-distance curves to get the corresponding mechanical properties.
- Tapping mode of AFM was adopted to scan the sample of graphene and Multi-wall Carbon Nanotube.
- The mechanical property of graphene—Young's Modulus, is calculated by in Hertz Model.



Methodology

Methods:

1. Graphene:

Ultrasound 20-30 min in the solution of propyl alcohol and water, to get a well-dispersion solution of graphene nanoplates.

2. Nanobubble:

Bubble was generated by electrolyze of electrolyte solution (KCl, 0.1 M) and then was scanned on the surface of button battery directly.

Principle of AFM:

The piezoelectric ceramic tube drives the AFM probe move in the vertical direction to maintain a constant interaction force between sample surface and AFM tip. The interaction force is detected by a laser reflected off the backside of the cantilever.

Results & Discussion

Graphene Nanoplate:




Fig. 1. The morphology of graphene nanoplate

Multi-wall Carbon Nanotube:




Fig. 2. The morphology of Carbon Nanotube by AFM

In this part, we perform the AFM characterization of Graphene Nanoplate and Carbon Nanotube, these two kind of nanomaterials have stacked-layer structure and tunnel-shape structure as shown in the figures respectively. The height results figure of Carbon Nanotube was plotted in MATLAB, as shown in Figure 2. (a).

Nanobubble:




Fig. 3. Approach-reside-retract cycle

- The bubbles were generated by electrical process, and AFM nanoindentation was performed under the liquid water.
- Hertz model was adopted to calculate the mechanical property of two spherical object. And the equation(left) can calculate the Young's Modulus.
- (Right) contact force-distance relationship of AFM's tip and nano-bubble. The MATLAB automatically deploy the 32*32 point of data, and integrate the data into the contour of nanobubble(bottom).

$$E = \frac{3(1 - \nu^2)P}{4} \left(\frac{1}{R_i} + \frac{1}{R} \right)^{\frac{1}{2}} h^{-\frac{3}{2}}$$





Fig. 4. (Up) 1024 contact line; (Bottom) Contour of nanobubble

Conclusion

- In this work, we get the basic information about Graphene Nanoplate and Multi-wall Carbon Nanotube, and those information includes the width distribution and the height of sample.
- Furthermore, the contact force and distance relationship has been obtained, and Young's Modulus could be calculated in the next step. More mechanical information about nanobubble would be reported in the future.

Acknowledgement:

The author thanks The University of Hong Kong (HKU) Summer Institute LABX 2019 Program.

Reference:

- Wiktorja Walczyk; Holger Schönherr; *Langmuir* **2014**, 30, 7112-7126.
- J. J. Roa1, G. O., J. Díaz, F. Sanz and M. Segarra, *Recent Patents on Nanotechnology* 2011, 5, 27-36.
- Kontomaris, S.-V., *Micro and Nanosystems*. 2018, 10, 11-22.